



READY-TO-USE DISINFECTANT CLEANER

- Kills 99.9% of bacteria in 30 seconds on hard nonporous surfaces
- Effective against Tuberculosis and Poliovirus
- Provides Hospital Disinfection. Meets the decontamination requirements of the OSHA Bloodborne Pathogens Standard
- Effective against HIV-1 (AIDS Virus), Rotavirus, Rhinovirus, RSV, Staph, Strep, Salmonella and Pseudomonas
- Citric acid formula contains no alcohol, bleach, phenol or harsh abrasives
- Will not dry out Naugahyde
- Proven compatibility with common surfaces in medical and dental offices, such as Naugahyde and vinyl chairs, laminate countertops, stainless steel sinks, plastic and porcelain
- No other product outperforms ¹
- For use in dental and physician offices, operating and emergency rooms, patient rooms, laboratories, nurses' stations, nursing homes, municipal facilities, schools, buses, ambulances and other vehicles
- No dilution necessary! Just pour into a clean, appropriately labeled trigger sprayer bottle
- EPA Reg. No. 675-55

¹ As determined by an independent testing facility*

Product	Order#	Size	Case Pack
Disinfectant Cleaner	36241-75691	1 gallon	4

Provides Efficacy Against the Following Microorganisms:

--

<p style="text-align: center;">Bacteria</p> <p>Salmonella choleraesuis Staphylococcus aureus Pseudomonas aeruginosa Enterococcus (Streptococcus) faecalis Escherichia coli Mycobacterium tuberculosis var bovis BCG Enterococcus faecium Streptococcus pyogenes Campylobacter jejuni</p>	<p style="text-align: center;">Fungi</p> <p>Trichophyton mentagrophytes</p> <p style="text-align: center;">Viruses</p> <p>Influenza A2 (Japan) Herpes Simplex Type 1 Herpes Simplex Type 2 Rhinovirus Type 39 Rotavirus Respiratory Syncytial Virus (RSV) Human Immunodeficiency Virus Type 1 (HIV-1) **(Aids Virus) Poliovirus Type 1 Echovirus 12 Adenovirus Type 2 Hepatitis A Virus Parainfluenza</p>
---	---